Standardisinç	g a sodiu	m hydroxide solution			
This is the example A2 of the EURACHEM / CITAC Guide "Quantifying Uncertainty in Analytical Measurement", Second Edition.					
A solution of sodium hydroxide (NaOH) is standardised against the titrimetric standard potassium hydrogen phtalate (KHP).					
The titrimetric standard KHP is dried and weighed. After the preparation of the NaOH solution a sample of the KHP is dissolved and then titrated using the NaOH solution.					
Model Equation {Molar mas					
M <sub>KHP</sub> = 8 *	M <sub>C</sub> + 5 * M <sub>t</sub>	<sub>H</sub> + 4 * M <sub>O</sub> + M <sub>K</sub> ;			
{Volume de	elivered by t	the piston burette, excluding repeatability}			
	-	on * f <sub>V-temperature</sub> ;			
(mass of K	HP weighed	d, excluding repeatability}			
$m_{KHP} = m_{cc}$	ontainer and KHF	p - M <sub>container less KHP</sub> ;			
{calculation of the concentration of the NaOH solution}					
c <sub>NaOH</sub> = (k <sub>mL</sub> * m <sub>KHP</sub> * P <sub>KHP</sub> ) / (M <sub>KHP</sub> * V <sub>T</sub> ) * f <sub>repeatability</sub> ;					
List of Quantiti					
Quantity	Unit	Definition			
M <sub>KHP</sub>	g/mol	Molar mass of KHP			
M <sub>C</sub>	g/mol	Atomic weight of carbon			
M <sub>H</sub>	g/mol	Atomic weight of hydrogen			
M <sub>o</sub>	g/mol	Atomic weight of oxygen			
M <sub>K</sub>	g/mol	Atomic weight of potassium			
V <sub>T</sub>	mL	Volume delivered by piston burette			
V <sub>nominal</sub>	mL	Volume indicated by burette			
f <sub>V</sub> -calibration		Uncertainty contribution to the volume due to uncertainty in calibratio the burette	n of		
f <sub>V-temperature</sub>		Uncertainty contribution to the volume due to temperature variation			
m <sub>KHP</sub>	g	Mass of KHP weighed			
m <sub>container and KHP</sub>	g	Mass of container and KHP			
m <sub>container less KHP</sub>					
С <sub>NaOH</sub>	mol/L				
k <sub>mL</sub>	mL/L	Conversion factor 1000 ml = 1L			
P <sub>KHP</sub>		Purity of the KHP			
f <sub>repeatability</sub>		Repeatability of the titration			
	Type B re				

Date: 08/29/2007	File: EURACHEM A2.smu	Page 1 of 4

	Standardising a sodium hydroxide solution				
The atomic weight of carbon and its uncertainty are taken from data listed in the latest IUPAC table of atomic weights. The IUPAC quoted data is considered to be of rectangular distribution.					
M <sub>H</sub> :	Type B rectangular distribution Value: 1.00794 g/mol Halfwidth of Limits: 0.00007 g/mol				
The atomic weight of hydrogen and its uncertainty are taken from data listed in the latest IUPAC table of atomic weights. The IUPAC quoted data is considered to be of rectangular distribution.					
M <sub>o</sub> :	Type B rectangular distribution Value: 15.9994 g/mol Halfwidth of Limits: 0.0003 g/mol				
	ght of oxigen and its uncertainty are taken from data listed in the latest IUF The IUPAC quoted data is considered to be of rectangular distribution.	PAC table of			
<b>М</b> <sub>К</sub> :	Type B rectangular distribution Value: 39.0983 g/mol Halfwidth of Limits: 0.0001 g/mol				
	ght of potassium and its uncertainty are taken from data listed in the latest The IUPAC quoted data is considered to be of rectangular distribution.	IUPAC table of			
V <sub>nominal</sub> :	Constant Value: 18.64 mL				
burette has three components, repeatability, calibration and temperature. The latter two are included in the uncertainty budget as separate factors. Repeatability of the volume delivery is taken into account via the combined repeatability term for the experiment, $f_{repeatability}$ . Another factor influencing the result of the titration, which can also be attributed to the automatic titration system, of which the burette is one part, is the bias of the end-point detection. The titration is performed under a protective atmosphere (Ar) to prevent absorption of CO <sub>2</sub> , which would bias the titration. There are no indications that the end-point determined from the shape of the pH-curve does not correspond to the equivalence-point, because a strong acid is titrated with a strong base. No further uncertainty contributions are introduced to cover the bias of the end-point detection.					
f <sub>V-calibration</sub> :	Type B triangular distribution Value: 1 Halfwidth of Limits: 0.0015				
The limits of accuracy for a 20 mL piston burette are indicated by the manufacturer as typically $\pm 0.03$ ml. No further statement is made about the level of confidence or the underlying distribution. An assumption is necessary to work with this uncertainty statement. In this case a triangular distribution is assumed. Since $f_{V-calibration}$ is a multiplicative factor to the nominal volume, which is only used to introduce the calibration uncertainty, it has the value 1. The halfwidth of limits corresponds to the relative uncertainty as stated by the manufacturer (i.e. 0.03 mL / 20 mL).					
f <sub>V-temperature</sub> :	Type B normal distribution Value: 1 Expanded Uncertainty: 0.0003 Coverage Factor: 1				
The laboratory temperature can vary by $\pm 3^{\circ}$ C. The uncertainty of the volume due to temperature variations can be calculated from the estimate of the possible temperature range and the coefficient of the volume expansion. The volume expansion of the liquid is considerably larger than that of the burette, so only the volume expansion of the liquid is considered. The coefficient of volume expansion for water is $2.1 \cdot 10^{-4}$ °C <sup>-1</sup> , which is used here also for the NaOH solution. This leads to a possible volume variation of $\pm (19 \cdot 3 \cdot 2.1 \cdot 10^{-4}/1.96)$ mL = $\pm 0.006$ mL. A rectangular distribution is assumed for the temperature					
Date: 08/29/2007	File: EURACHEM A2.smu	Page 2 of 4			

<b>n</b> <sub>container and</sub>	<sub>КНР</sub> : Type B rectang Value: 60.5450 Halfwidth of Lin	g		
systematic o contributing quotes ±0.1 convert this	offset across the scale source of uncertainty is 5 mg for the linearity. T linearity contribution in	will also cancel due to t s the linearity of the bal he manufacturer recom to a standard uncertatir	e combined repeatability term, f <sub>repeatability</sub> . Any he wheighing by difference. The only ance. The calibration certficate of the balance mends using a rectangular distribution to ny. This uncertainty is not correlated with the observations and the linearity effects are not	
<b>n</b> <sub>container less</sub>	<sub>КНР</sub> : Туре В rectang Value: 60.1562 Halfwidth of Lin	g		
systematic o contributing quotes ±0.1 convert this	offset across the scale source of uncertainty is 5 mg for the linearity. T linearity contribution in	will also cancel due to t s the linearity of the bal he manufacturer recom to a standard uncertatir	e combined repeatability term, f <sub>repeatability</sub> . Any he wheighing by difference. The only ance. The calibration certficate of the balance mends using a rectangular distribution to ny. This uncertainty is not correlated with the observations and the linearity effects are not	
« <sub>mL</sub> :	Constant Value: 1000 mL	/L		
P <sub>KHP</sub> : Type B rectangular distribution Value: 1 Halfwidth of Limits: 0.0005				
			as $1.0000 \pm 0.0005$ . No further information assumed to be of rectangular distribution.	
repeatability	Type B normal Value: 1 Expanded Unce Coverage Facto	ertainty: 0.0005		
ncludes at l nagnitude c shows that t	east the repeatability o of this uncertainty contr he overall repeatability result, which is only us of 0.0005.	f the wheighings and of ibution is assessed dur of the titration experim	the operations are combined in this factor. It the volume delivered by the burette. The ng the method validation stage. The data ent is 0.05%. Since $f_{repeatability}$ is a multiplicative eatability uncertainty, it has the value 1 with an	
Quantity	Value	Standard Uncertainty		
M <sub>KHP</sub>	204.221200 g/mol	3.765.10 <sup>-3</sup> g/mol		
V <sub>T</sub>	18.64000 mL	0.01271 mL	1	
m <sub>KHP</sub>	0.3888000 g	122.5·10 <sup>−6</sup> g	]	
			Page 3 of 4	

Standardising a sodium hydroxide solution

NaOH	Concentration of	f the sodium h	ydroxide sol	ution		
Quantity	Value	Standard Uncertainty	Distributio n	Sensitivity Coefficient	Uncertainty Contribution	Index
M <sub>KHP</sub>	204.221200 g/mol	3.765∙10 <sup>-3</sup> g/mol				
M <sub>C</sub>	12.0107000 g/mol	461.9∙10 <sup>-6</sup> g/mol	rectangular	-4.0·10 <sup>-3</sup>	-1.8·10 <sup>-6</sup> mol/L	0.0 %
M <sub>H</sub>	1.00794000 g/mol	40.41∙10 <sup>-6</sup> g/mol	rectangular	-2.5·10 <sup>-3</sup>	-100·10 <sup>-9</sup> mol/L	0.0 %
M <sub>O</sub>	15.9994000 g/mol	173.2·10 <sup>-6</sup> g/mol	rectangular	-2.0·10 <sup>-3</sup>	-350·10 <sup>-9</sup> mol/L	0.0 %
M <sub>K</sub>	39.09830000 g/mol	57.74·10 <sup>-6</sup> g/mol	rectangular	-500·10 <sup>-6</sup>	-29·10 <sup>-9</sup> mol/L	0.0 %
V <sub>T</sub>	18.64000 mL	0.01271 mL				
V <sub>nominal</sub>	18.64 mL					
f <sub>V-calibration</sub>	1.0000000	612.4·10 <sup>-6</sup>	triangular	-0.10	-63·10 <sup>-6</sup> mol/L	41.8 %
f <sub>V-temperature</sub>	1.0000000	300.0·10 <sup>-6</sup>	normal	-0.10	-31.10 <sup>-6</sup> mol/L	10.0 %
m <sub>KHP</sub>	0.3888000 g	122.5∙10 <sup>-6</sup> g				
n <sub>container</sub> and KHP	60.54500000 g	86.60·10 <sup>-6</sup> g	rectangular	0.26	23.10 <sup>-6</sup> mol/L	5.5 %
n <sub>container less KHP</sub>	60.15620000 g	86.60·10 <sup>-6</sup> g	rectangular	-0.26	-23·10 <sup>-6</sup> mol/L	5.5 %
k <sub>mL</sub>	1000.0 mL/L					
P <sub>KHP</sub>	1.0000000	288.7 <b>·</b> 10 <sup>-6</sup>	rectangular	0.10	29-10 <sup>-6</sup> mol/L	9.3 %
f <sub>repeatability</sub>	1.0000000	500.0·10 <sup>-6</sup>	normal	0.10	51.10 <sup>-6</sup> mol/L	27.8 %
C <sub>NaOH</sub>	0.10213616 mol/L	96.78·10 <sup>-6</sup> mol/L				

## **Results:**

Quantity	Value	Expanded Uncertainty	Coverage factor	Coverage
с <sub>NaOH</sub>	0.10214 mol/L	190∙10 <sup>-6</sup> mol/L	2.00	95% (t-table 95.45%)